


For Research Use Only (RUO)
Not for use in diagnostic procedures

● CHROMagar™ **Acinetobacter**



For detection of *Acinetobacter*
and MDR *Acinetobacter* species

CHR  **Magar**
The Chromogenic Media Pioneer

● CHROMagar™ Acinetobacter

For Research Use Only (RUO). Not for use in diagnostic procedures.





Plate Reading

For detection of *Acinetobacter* sp.:

- *Acinetobacter* sp.
→ red
- Other gram (-)
→ mostly inhibited
- Gram(+) bacteria and yeasts
→ inhibited

For detection of MDR *Acinetobacter* sp.
(if using the optional supplement CR102):

- MDR *Acinetobacter*
→ red
- Non-MDR *Acinetobacter*
→ inhibited
- Other gram (-)
→ mostly inhibited
- Gram(+) bacteria and yeasts
→ inhibited



Quality Control Strains

Acinetobacter baumannii ATCC® 19606 red
Acinetobacter calcoaceticus ATCC® 51432 red
E. faecalis ATCC® 29212 inhibited
C. tropicalis ATCC® 1369 inhibited

using CR102 supplement:

Acinetobacter baumannii ATCC® 19606 inhibited

ATCC® is a registered trademark of the American Type Culture Collection

For detection of *Acinetobacter* and MDR *Acinetobacter* sp.

Medium Description

Powder Base	Total	32.8 g/L
	Agar	15.0
+	Peptone and yeast extract	12.0
	Salts	4.0
+	Chromogenic mix	1.8
	Storage at 15/30°C - pH: 7.0 +/-0.2	
	Shelf Life	2 years
Supplement (included in the pack)	Growth and regulator factors	4ml/L
	Storage at 15/30°C	
+	Aspect: Liquid Form	
	Shelf Life	2 years
CHROMagar MDR Supplement : CR102 Order separately	Selective mix	5 doses
	(1 dose qsf 1000ml of final media)	
	Storage at 2/8°C	
	Shelf Life	2 years

Please read carefully the instructions for use (IFU document) available on www.drg-international.com.

For US Customers: for Research Use Only

Order References

Please use these product references when contacting DRG

5000 ml pack AC092 (Product = AC092(B) powder base + AC092(S) supplement)

CHROMagar MDR (Optional) : 5000 ml pack.....CR102

Manufacturer: CHROMagar
4 place du 18 juin 1940 75006 Paris - France

Distributed in the US by DRG International
www.DRG-International.com