

# CHROMagar™ **Y.enterocolitica**

**For detection and direct differentiation  
of pathogenic *Y.enterocolitica***

# ● CHROMagar™ *Y. enterocolitica*

www.CHROMagar.com

## Plate Reading

For detection of:

- Pathogenic *Y. enterocolitica*  
→ Mauve
- Non pathogenic *Y. enterocolitica* and background flora (*Citrobacter*, *Enterobacter*, *Aeromonas* etc)  
→ Metallic blue, inhibited or limited growth in metallic blue colour

## For detection and direct differentiation of pathogenic *Yersinia enterocolitica* strains

### Background

Among the *Yersinia* genus, *Yersinia enterocolitica* is one of the most common food borne pathogen. In several countries, *Y. enterocolitica* has eclipsed *Shigella* and approaches *Salmonella* and *Campylobacter* as the predominant cause of acute bacterial gastroenteritis. Its ability to grow at refrigeration temperature makes it an increasing concern in terms of food safety. This germ most commonly affects young individuals. However, only a few strains of *Y. enterocolitica* cause illness in humans. Those pathogenic *Y. enterocolitica* strains belong to biotypes 1B, 2, 3, 4, and 5, whereas biotype 1A strains are non-pathogenic and widespread in the environment. The major animal *Y. enterocolitica* reservoir causing illnesses are pigs.

### Medium Performance

#### 1 HIGH SPECIFICITY OF THE MAUVE COLOUR

##### 1 Clear differentiation of *Yersinia* among background flora


Background flora, such as *Citrobacter*, often has a *Yersinia*-like red aspect on CIN agar which generate a heavy routine workload on false positive results. Thanks to CHROMagar *Y. enterocolitica*, this unnecessary workload will be eliminated as only pathogenic *Yersinia enterocolitica* colonies will grow mauve.

##### 2 Colourful differentiation of pathogenic among non pathogenic *Yersinia enterocolitica*

In traditional culture media like CIN agar, pathogenic and non-pathogenic biotypes have the same aspect. On CHROMagar *Y. enterocolitica*, the pathogenic biotypes grow in a distinctive mauve colour, differentiated at a glance, from the non pathogenic which will grow metallic blue.

#### 2 HIGH SELECTIVITY: REDUCTION OF BACKGROUND FLORA

Background flora is dramatically reduced on CHROMagar *Y. enterocolitica* allowing an easy reading of the plates. **The laboratory will concentrate its efforts and resources only on suspect colonies that have a real potential of pathogenicity.**

### Medium Description

<b>Powder Base</b>	Total ..... 41.3 g/L Agar ..... 15.0 Peptones ..... 20.0 Salts ..... 5.0 Chromogenic & selective mix ..... 1.3 Storage at 15/30°C - pH: 7.0 ± 0.2 Shelf Life ..... 2 years
<b>+</b> <b>Supplement</b> (Included in the pack)	Powder form ..... 100 mg/L Storage at 2/8°C      Shelf Life ..... 18 months

Usual Samples	All
Procedure	Direct Streaking or after appropriate enrichment step. Incubation at 30°C, 24h-36h. Aerobic conditions.
Please read carefully the instructions for use (IFU document) available on <a href="http://www.CHROMagar.com">www.CHROMagar.com</a>	

## Quality Control Strains

<i>Y. enterocolitica</i> pYV+ ATCC® 23715	..... mauve
<i>Y. enterocolitica</i> pYV-	..... metallic blue
<i>E. coli</i> ATCC® 25922	..... inhibited
<i>E. faecalis</i> ATCC® 29212	..... inhibited
<i>Pseudomonas</i> ATCC® 9027	..... inhibited
<i>Citrobacter</i> ATCC® 8090	..... inhibited

ATCC® is a registered trademark of the American Type Culture Collection


## ORDER REFERENCES

Please use these product references when contacting your local distributor:

5000 ml pack ..... YE492  
(Included in that reference: powder base YE492(B) + supplement YE492(S) )

Manufacturer: CHROMagar  
4 place du 18 juin 1940 75006 Paris - France  
Email: [CHROMagar@CHROMagar.com](mailto:CHROMagar@CHROMagar.com)  
Website: [www.CHROMagar.com](http://www.CHROMagar.com)  
Find your nearest distributor on  
[www.CHROMagar.com/contact](http://www.CHROMagar.com/contact)

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