CHROMagar™ **Listeria**

Instructions For Use

Available in several languages

NT-EXT-009

Version 5

ENGLISH

English Version

FRANCAIS

Version Française

ESPAÑOL


Version Español

DEUTSCH

Deutsch Version

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MEDIUM PURPOSE

Chromogenic medium for detection, isolation and enumeration of *L.monocytogenes*

Listeria monocytogenes is a widespread bacteria, present in the soil, sewage or faecal matter. Its ability to form listerial biofilms on contact surfaces makes it difficult to eliminate. This pathogen can cause serious food poisoning and is therefore frequently a microbial Q.C. target in food processing facilities to avoid food contamination. Contamination can occur at all steps of the food manufacturing chain from raw materials to place of consumption.

COMPOSITION

The product is composed of a powder base (B) and a supplement (S).

Product	=	Base (B)	+	Supplement (S)
Total g/L		51.5 g/L		9.0 g/L
Composition g/L		Agar 15.0 Peptone and yeast extract 23.0 NaCl 5.0 Chromogenic mix 8.5		Selective and enrichment mix 9.0
Aspect		Powder Form		Powder Form
STORAGE		2/30°C		2/8°C
FINAL MEDIA pH		7.0 +/- 0.5		

APPLICATIONS

• In the food field:

For use with all human food products and environmental samples:

CHROMagar Listeria method (presence / absence):

→ Certified NF VALIDATION CHR-21/1-12/01.

Validity Date: 14-Dec-2017

CHROMagar Listeria enumeration method:

→ Certified NF VALIDATION CHR-21/2-12/06.

Validity Date: 14-Dec-2017

• In the clinical field:

For clinical field, allowing isolation of *L.monocytogenes* in samples of human origin eg. blood cultures, vaginal swabs, or CSF.

PREPARATION (Calculation for 1L)

Step 1

Preparation of the base
CHROMagar
Listeria base (B)

- Disperse slowly 51.5g of powder base in 1L of purified water.
- Stir until agar is well thickened.
- Heat to 121°C +/- 1°C for 15 minutes.
- Cool in a water bath to 47°C +/- 2°C.

Step 2

Preparation of the
CHROMagar
Listeria supplement (S)


- Add 9g to 40ml of STERILE purified water.
- Add a magnetic bar and stir at high speed (700-1000 rpm) at least 30 minutes WITHOUT HEATING until getting a creamy homogeneous suspension as below:

Final Media HELPING CALCULATION

1 L --> Add 9g of supplement in 40 ml of sterile water

5 L --> Add 45g of supplement in 200 ml of sterile water

STAGE 1:
Possible Presence of grains on the surface of the vial and in the mix

Aspect at this step:
• Liquid with undissolved grains

STAGE 2: Formation of white dense foam on the surface and rise of foam

Aspect at this step:
• Dense/Compact foam (SURFACE)
• Foamy Mix

STAGE 3: Less dense/compact Foam, more creamy mix

Aspect at this step:
• No Rise of Foam
• More ventilated foam
• Liquid

STAGE 4:
VERY GOOD MIX

Aspect at this step:
• Homogeneous
• Fluid (liquid)
• Without grains No Foam
GOOD PREPARATION

Step 3

Mixing of the prepared base (B) and the prepared supplement (S)

- Place the melted 47°C cooled CHROMagar Listeria base under gentle stirring.
- Add the homogeneous reconstituted supplement, keeping the gentle stirring during 1 or 2 minutes until complete homogenisation.
- Pour IMMEDIATELY into sterile Petri dishes.

Warning: DO NOT STACK THE PETRI DISHES

- Let them cool down. Let dry.

Storage

- Store in the dark before use.
- Prepared media plates can be kept for one day at room temperature.
- Plates can be stored for up to two weeks under refrigeration (2/8°C) if properly prepared and protected from light and dehydration.

CHROMagar™ Listeria

INOCULATION

Please refer to Part **Procedures**

- If the agar plate has been refrigerated, allow to warm to room temperature before inoculation.

Collection and transport of the samples should be well handled and adapted to the particular specimen according to good laboratory practices.

PROCEDURES

CHROMagar Listeria Method

• **For the Detection of Presence / Absence of *L.monocytogenes* in all human food products and environmental samples. (Refer to Illustration 1)**

- 1 Perform an enrichment step of the sample, mixing \underline{x} g or \underline{x} ml in $9\underline{x}$ Fraser 1/2 broth
→ 24h +/- 2h at 30°C +/- 1°C
- 2 Spread 0,1ml on a well dried CHROMagar Listeria plate, by streaking on 1/3 of the plate until liquid is adsorbed. Isolate on the rest of the plate and incubate.
→ 24h +/- 2h at 37°C +/- 1°C
- 3 Typical colonies:
NO → Absence of *L.monocytogenes*
YES → In the context of NF Validation, all samples identified as positive by the alternative method must be confirmed in one of the following ways:
- 4 1/ using the conventional tests described in the standardized methods by CEN or ISO (including the purification step).
2/ with the CHROMagar Identification Listeria directly from the suspect colony without step purification.

In the event of discordant results (positive presumptive by the alternative method, non-confirmed by one of the means described above) the laboratory must follow the necessary steps to ensure validity of the result obtained.

Test portions weighing more than 25 g have not been tested.

For environmental samples, typical colonies should be absent after 24h, incubate for additional 24h +/- 2h.

CHROMagar Listeria enumeration Method

• **For Enumeration of *L.monocytogenes* in all human food products and environmental samples. (Refer to Illustration 2)**

- 1 Perform a revivication step to the sample and mix \underline{x} g or \underline{x} ml in $9\underline{x}$ BPW or non selective half Fraser broth
→ 1h +/- 5 min at 20°C +/- 1°C
- 2 Spread 0,1 ml on one well dried CHROMagar Listeria plate (in case of estimation of small numbers spread 1ml on 3 plates) and incubate.
→ 48h +/- 2h at 37°C +/- 1°C
A first reading at 24h allows a more rapid detection of heavily contaminated samples. However the final result count is reached after 48 h +/- 2h. Refer to ISO 7218 standard for calculation and interpretation of results. The spread of a single dilution may present a risk of gap counts compared to the reference method.
- 3 Typical colonies:
NO → Absence of *L.monocytogenes*
YES → Enumerate and confirm a colony directly from CHROMagar Listeria. If you do not confirm 5 colonies counting, it can create a risk of making a overestimated result because of the possible presence of typical colonies that would not be *L. monocytogenes*. In the context of NF Validation: If confirmation has already been done during the detection phase with the CHROMagar Listeria method, a new confirmation is not necessary after the CHROMagar Listeria enumeration method. Otherwise, all samples identified as positive by the alternative method must be confirmed in one of the following ways:
- 4 1/ using the conventional tests described in the standardized methods by CEN or ISO (including the purification step).
2/ with the CHROMagar Identification Listeria directly from the suspect colony without step purification.

In the event of discordant results (positive presumptive by the alternative method, non-confirmed by one of the means described above) the laboratory must follow the necessary steps to ensure validity of the result obtained.

Illustration 2: Enumeration Method



Illustration 1: Detection Method


CHROMagar™ Listeria

INTERPRETATION

Microorganism	Typical colony appearance
<i>L.monocytogenes</i>	→ blue, diameter less than 3mm, regular and white halo
Other microorganisms	→ blue, colourless, other colour, inhibited

Typical colony appearance on CHROMagar Listeria

The pictures shown are not contractual.

PERFORMANCE & LIMITATIONS

• Some strains of *L.ivanovii* may also give blue colonies with white halo distinguished during identification. The species *L.ivanovii* is rarely found in food. Some strains of *B.cereus* can also grow as blue colonies. They can easily be distinguished from colonies of *L.monocytogenes* as they are much larger with an irregular edge to the colony and very large white halo.

QUALITY CONTROL

Please perform Quality Control according to the use of the medium and the local QC regulations and norms. Good preparation of the medium can be tested, isolating the ATCC strains below:

Microorganism	Typical colony appearance
<i>L. monocytogenes</i> ATCC® 19115*	→ blue with halo
<i>L. innocua</i> CIP 8012 = ATCC® 33091*	→ blue without halo
<i>E. faecalis</i> ATCC® 29212 = WDCM 00087	→ inhibited
<i>E. coli</i> ATCC® 25922 = WDCM 00013	→ inhibited

WDCM : World Data Center for Microorganism World Federation for Culture Collection

* Corresponding strains on www.straininfo.net

WARNINGS

- Do not use plates if they show any evidence of contamination or any sign of deterioration.
- Do not use the product beyond its expiry date or if product shows any evidence of contamination or any sign of deterioration.
- For *in vitro* diagnostic use. This laboratory product should be used only by trained personnel in compliance with good laboratory practices.
- Any change or modification in the procedure may affect the results.
- Any change or modification of the required storage temperature may affect the performance of the product.
- Unappropriate storage may affect the shelf life of the product.
- Recap the bottles/vials tightly after each preparation and keep them in a low humidity environment, protected from moisture and light.
- For a good microbial detection: collection and transport of specimen should be well handled and adapted to the particular specimen according to good laboratory practices.

DISPOSAL OF WASTE

After use, all plates and any other contaminated materials must be sterilized or disposed of by appropriate internal procedures and in accordance with local legislations. Plates can be destroyed by autoclaving at 121°C for at least 20 minutes.

REFERENCES

Please refer to our website page «Publications» for scientific publications about this particular product.

Web link: <http://www.chromagar.com/publication.php>

IFU/LABEL INDEX

- Quantity of powder sufficient for X liters of media
- Expiry date
- Required storage temperature
- Store away from humidity

Need some Technical Documents?

Available for download on www.CHROMagar.com

- Certificate of Analysis (CoA) --> One per Lot
- Material Safety Data Sheet (MSDS)

Pack Size

Ordering References

Base (B)

Supplement (S)

1000 ml

50 Tests of 20ml

=

LM851

=

LM851(B)
Weight: 51.5gr

+

LM851(S)
Weight: 9gr

5000 ml

250 Tests of 20ml

=

LM852

=

LM852(B)
Weight: 257.5gr

+

LM852(S)
Weight: 45gr

NT-EXT-009 V5 / 18-Oct-13

CHROMagar is a trademark created by Dr A. Rambach

ATCC® is a registered trademark of the American Type Culture Collection .

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